




USER MANUAL – EN

IN 23198 Facial steamer inSPORTline Cison






CONTENTS

SAFETY INSTRUCTIONS..... 3
PRODUCT DESCRIPTION 3
ASSEMBLY 4
USE 4
STEPS OF USE 5
CLEANING AND MAINTENANCE 5
SPECIFICATIONS..... 5
ENVIRONMENT PROTECTION 6
TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS..... 6

SAFETY INSTRUCTIONS

- Read the manual before use and keep it for future reference.
- Always follow all instructions, warnings, and cautions.
- Do not use the product near water or in places with high humidity.
- Not suitable for people sensitive to heat.
- If you experience pain or observe any skin abnormalities, stop using the product immediately.
- The device is not suitable for: pregnant women, people with sensitive skin, suffering from allergies, burns, dermatological diseases, people recovering from surgery, etc.
- Do not touch the metal parts of the plug when connecting.
- Switch off and unplug after use.
- Do not connect or disconnect the product with wet hands.
- The product is not a toy, keep out of reach of children and pets.
- Do not use near flammable materials (sprays, gas, etc.)
- Allow to cool after use and then pour out the rest of the water.
- Use only distilled or filtered water.
- During use, the device emits hot steam, keep the appliance at least 20 cm away from your body.
- Do not shake, tilt, or carry the device during use or for 10 minutes after use. Hot water could spill and cause scalding.
- During using, place the product on a flat surface, do not hold it in your hands.
- Do not use more than twice a day for 10 minutes. For normal skin, use 3-5 times a week, for sensitive skin once a week.
- If you experience any discomfort, stop using the product immediately.

PRODUCT DESCRIPTION


	1.	Steamer
	2.	Switch
	3.	The body of the device

ASSEMBLY

Place the device on a flat and solid surface.

Measure out the water, do not exceed the marked gauge.

Align the steamer with the main body of the device, there is a mark (Pic. 2), turn it 90 ° clockwise. The device is ready for use (Pic. 3).



USE


Use only distilled water or filtered water.

Due to the different quality of tap water, we recommend using only distilled or filtered water.

Do not mix with other liquids, oils, etc.

1. Plug in the device.
2. The steam starts to evaporate 20 seconds after turning on. If you overfill the container, there is a risk that water droplets will escape with the steam. Do not exceed the marked gauge.
3. The steamer must be at least 20 cm away from the face. If the tank starts to run out of water, the steam will be weaker. As soon as the water is depleted, the device switches off. If you want to continue using the device, turn it off, wait 10 minutes, and repeat steps 1 - 3. Otherwise, turn off the device, wait for it to cool, and store it.
4. Turn off the device after use. Then perform a routine skin care.

STEPS OF USE



1. Pour distilled / filtered water
2. Attach the head
3. Turn on the device

CLEANING AND MAINTENANCE

After use, pour out all the water and wait for the device to dry.

1. Allow the product to cool down for at least 10 min.
2. Turn off the device, pour off the excess water and dry the surface of the device.
3. When draining the device, unplug the steamer and pour out the water from the body of the device.

Clean the product from dust at least once a month.

Add water to the water tank, attach the steamer and shake gently 3x - 4x, do not tilt the product while shaking. Then pour the water out of the tank. Repeat 2x - 3x.

SPECIFICATIONS

Name:	Nano Ionic Facial Steamer
Rated power:	300 W
Rated voltage:	110V – 220V
Steam volume:	6.5g / min
Tank capacity:	65 ml

ENVIRONMENT PROTECTION

After the product lifespan expired or if the possible repairing is uneconomic, dispose it according to the local laws and environmentally friendly in the nearest scrapyards.

By proper disposal you will protect the environment and natural sources. Moreover, you can help protect human health. If you are not sure in correct disposing, ask local authorities to avoid law violation or sanctions.

Don't put the batteries among house waste but hand them in to the recycling place.

TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS

General Conditions of Warranty and Definition of Terms

All Warranty Conditions stated here under determine Warranty Coverage and Warranty Claim Procedure. Conditions of Warranty and Warranty Claims are governed by Act No. 89/2012 Coll. Civil Code, and Act No. 634/1992 Coll., Consumer Protection, as amended, also in cases that are not specified by these Warranty rules.

The seller is SEVEN SPORT s.r.o. with its registered office in Strakonická street 1151/2c, Prague 150 00, Company Registration Number: 26847264, registered in the Trade Register at Regional Court in Prague, Section C, Insert No. 116888.

According to valid legal regulations it depends whether the Buyer is the End Customer or not.

"The Buyer who is the End Customer" or simply the "End Customer" is the legal entity that does not conclude and execute the Contract in order to run or promote his own trade or business activities.

"The Buyer who is not the End Customer" is a Businessman that buys Goods or uses services for the purpose of using the Goods or services for his own business activities. The Buyer conforms to the General Purchase Agreement and business conditions.

These Conditions of Warranty and Warranty Claims are an integral part of every Purchase Agreement made between the Seller and the Buyer. All Warranty Conditions are valid and binding, unless otherwise specified in the Purchase Agreement, in the Amendment to this Contract or in another written agreement.

Warranty Conditions

Warranty Period

The Seller provides the Buyer a 24 months Warranty for Goods Quality, unless otherwise specified in the Certificate of Warranty, Invoice, Bill of Delivery or other documents related to the Goods. The legal warranty period provided to the Consumer is not affected.

By the Warranty for Goods Quality, the Seller guarantees that the delivered Goods shall be, for a certain period of time, suitable for regular or contracted use, and that the Goods shall maintain its regular or contracted features.

The Warranty does not cover defects resulting from (if applicable):

- User's fault, i.e. product damage caused by unqualified repair work, improper assembly, insufficient insertion of seat post into frame, insufficient tightening of pedals and cranks
- Improper maintenance
- Mechanical damages
- Regular use (e.g. wearing out of rubber and plastic parts, moving mechanisms, joints, wear of brake pads/blocks, chain, tires, cassette/multi wheel etc.)
- Unavoidable event, natural disaster
- Adjustments made by unqualified person
- Improper maintenance, improper placement, damages caused by low or high temperature, water, inappropriate pressure, shocks, intentional changes in design or construction etc.

Warranty Claim Procedure

The Buyer is obliged to check the Goods delivered by the Seller immediately after taking the responsibility for the Goods and its damages, i.e. immediately after its delivery. The Buyer must check the Goods so that he discovers all the defects that can be discovered by such check.

When making a Warranty Claim the Buyer is obliged, on request of the Seller, to prove the purchase and validity of the claim by the Invoice or Bill of Delivery that includes the product's serial number, or eventually by the documents without the serial number. If the Buyer does not prove the validity of the Warranty Claim by these documents, the Seller has the right to reject the Warranty Claim.

If the Buyer gives notice of a defect that is not covered by the Warranty (e.g. in the case that the Warranty Conditions were not fulfilled or in the case of reporting the defect by mistake etc.), the Seller is eligible to require a compensation for all the costs arising from the repair. The cost shall be calculated according to the valid price list of services and transport costs.

If the Seller finds out (by testing) that the product is not damaged, the Warranty Claim is not accepted. The Seller reserves the right to claim a compensation for costs arising from the false Warranty Claim.

In case the Buyer makes a claim about the Goods that is legally covered by the Warranty provided by the Seller, the Seller shall fix the reported defects by means of repair or by the exchange of the damaged part or product for a new one. Based on the agreement of the Buyer, the Seller has the right to exchange the defected Goods for a fully compatible Goods of the same or better technical characteristics. The Seller is entitled to choose the form of the Warranty Claim Procedures described in this paragraph.

The Seller shall settle the Warranty Claim within 30 days after the delivery of the defective Goods, unless a longer period has been agreed upon. The day when the repaired or exchanged Goods is handed over to the Buyer is considered to be the day of the Warranty Claim settlement. When the Seller is not able to settle the Warranty Claim within the agreed period due to the specific nature of the Goods defect, he and the Buyer shall make an agreement about an alternative solution. In case such agreement is not made, the Seller is obliged to provide the Buyer with a financial compensation in the form of a refund.

CZ

SEVEN SPORT s.r.o.

Registered Office: Strakonická 1151/2c, Praha 5, 150 00, ČR
Headquarters: Dělnická 957, Vítkov, 749 01
Warranty & Service: Čermenská 486, Vítkov 749 01

CRN: 26847264
VAT ID: CZ26847264

Phone: +420 556 300 970
E-mail: eshop@insportline.cz
reklamace@insportline.cz
servis@insportline.cz

Web: www.inSPORTline.cz

SK

inSPORTline s.r.o.

Headquarters, warranty & service center: Električná 6471, Trenčín 911 01, SK

CRN: 36311723
VAT ID: SK2020177082

Phone: +421(0)326 526 701
E-mail: objednavky@insportline.sk
reklamacie@insportline.sk
servis@insportline.sk

Web: www.inSPORTline.sk