



USER MANUAL – EN
IN 26060 Cupping set inSPORTline Sebsahar



Dimensions:

2 ks, $\Phi 7.5 \times \Phi 5.3 \times 9.6$ cm

4 ks, $\Phi 6.5 \times \Phi 4.5 \times 8.5$ cm

4 ks, $\Phi 5.8 \times \Phi 4.0 \times 8.0$ cm

2 ks, $\Phi 5.0 \times \Phi 3.2 \times 7.2$ cm

2 ks, $\Phi 5.0 \times \Phi 3.2 \times 6.3$ cm

SevenSport s.r.o. reserves the right to make any changes and improvements to its product without prior notice. Visit our website www.insportline.eu where you will find the latest version of the manual.

SAFETY INSTRUCTIONS


- Read the manual before use and keep it for future reference.
- The manufacturer is not responsible for damage or injury to health caused by improper use.
- Always perform cupping after consulting a doctor.
- Do not apply the cups to the areas of: the heart and large blood vessels, ears, eyes, nose, nipples, genitals, varicose veins, damaged skin, areas of pigment spots, for pregnant women the abdomen and loins
- Not suitable for people suffering from:
 - Inflammation or skin irritation
 - Increased fragility of the veins
 - Any disease or inflammation of blood vessels
 - With a blood clotting disorder (e.g. hemophilia)
 - Shortness of breath
 - Autoimmune disease
 - Multiple sclerosis
 - Anemia and general cachexia
 - High fevers (above 38.5 °C)
 - Unstable hypertension, heart disease
 - Epilepsy
 - Viral disease
 - After surgery and transplantation
 - Pregnant women may be cupped only after 4 months (there is a risk of miscarriage earlier), and from the 5th month of pregnancy, the cups may not be applied to the abdomen and loins. Strong pressure must not be used during pregnancy
 - Tumor disease
 - After consuming alcohol or drugs
 - Do not apply immediately after a meal

WARNING: Patients taking medications related to anticoagulant therapy or blood pressure must consult a physician for use.

WARNING: The patient should feel warmth and stretching, not pain. If the patient feels pain, remove the cup immediately.

WARNING: In some cases (if the cup is suctioned for a long time or too strongly), blisters may appear on the skin. In this case, remove the flask immediately.

USE

	1. neck and feet stimulation
	2. stimulation of the neck, feet, hands, legs
	3. stimulation of shoulders, abdomen, legs and back
	4. stimulation of the shoulders, abdomen, back, waist and thighs
	5. stimulation

BENEFITS



USE

PREPARATION

Depending on the place of use, select two suitable cups (one for replacement). Prepare a cup, tweezers, a container of 95% alcohol (with a large opening for dipping cotton wool), cotton wool, matches, a towel, soap and a sink.

CHECK BEFORE APPLICATION

According to the established diagnosis, check the correct placement of the cup. Check that there are no sharp edges on the cup.

CUP APPLICATION

Clean the skin. Grab a cotton ball with tweezers and soak it in 95% alcohol. Light the soaked cotton ball and hold it in the flask for 3 seconds, then pull the cotton ball out with a quick movement and immediately press the cup against the skin. The flame sucks the air out of the cup and creates a vacuum that allows the cup to be suctioned to the skin.

WARNING: The patient should feel warmth and stretching, not pain. If the patient feels pain, remove the cup immediately.

WARNING: In some cases (if the cup is suctioned for a long time or too strongly), blisters may appear on the skin. In this case, remove the flask immediately.

LENGTH OF APPLICATION

Apply for 10-20 min (for healthy adult person)

CUP REMOVING

If the cup is released on its own and the skin does not show any cracking or damage to the small capillaries, it is possible to apply the cup again.

To remove, gently squeeze the skin at the edge of the cup. The cup takes in air, releases painlessly and can be removed. The patient should then lie covered and rest in the warm environment for 15-20 minutes.

AMOUNT OF ALCOHOL

More alcohol will cause more vacuum, less alcohol less vacuum. Always adjust the amount as needed.

USE INTERVAL

It is advisable to consult a professional about the frequency and duration of use, the time varies fundamentally depending on the type of use.

LOCATIONS

Shoulders, chest, back, waist, hips, ribs, cervical vertebrae, ankles, calves and other muscles. In addition, you can use the flasks on other parts of the body, as needed.

ENVIRONMENT PROTECTION

After the product lifespan expired or if the possible repairing is uneconomic, dispose it according to the local laws and environmentally friendly in the nearest scrapyard.

By proper disposal you will protect the environment and natural sources. Moreover, you can help protect human health. If you are not sure in correct disposing, ask local authorities to avoid law violation or sanctions.

Don't put the batteries among house waste but hand them in to the recycling place.

TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS

General Conditions of Warranty and Definition of Terms

All Warranty Conditions stated here under determine Warranty Coverage and Warranty Claim Procedure. Conditions of Warranty and Warranty Claims are governed by Act No. 89/2012 Coll. Civil Code, and Act No. 634/1992 Coll., Consumer Protection, as amended, also in cases that are not specified by these Warranty rules.

The seller is SEVEN SPORT s.r.o. with its registered office in Strakonická street 1151/2c, Prague 150 00, Company Registration Number: 26847264, registered in the Trade Register at Regional Court in Prague, Section C, Insert No. 116888.

According to valid legal regulations it depends whether the Buyer is the End Customer or not.

“The Buyer who is the End Customer” or simply the “End Customer” is the legal entity that does not conclude and execute the Contract in order to run or promote his own trade or business activities.

“The Buyer who is not the End Customer” is a Businessman that buys Goods or uses services for the purpose of using the Goods or services for his own business activities. The Buyer conforms to the General Purchase Agreement and business conditions.

These Conditions of Warranty and Warranty Claims are an integral part of every Purchase Agreement made between the Seller and the Buyer. All Warranty Conditions are valid and binding, unless otherwise specified in the Purchase Agreement, in the Amendment to this Contract or in another written agreement.

Warranty Conditions

Warranty Period

The Seller provides the Buyer a 24 months Warranty for Goods Quality, unless otherwise specified in the Certificate of Warranty, Invoice, Bill of Delivery or other documents related to the Goods. The legal warranty period provided to the Consumer is not affected.

By the Warranty for Goods Quality, the Seller guarantees that the delivered Goods shall be, for a certain period of time, suitable for regular or contracted use, and that the Goods shall maintain its regular or contracted features.

Batteries

6-month battery warranty – we guarantee that battery's nominal capacity does not fall below 70% of its total capacity within 6 months of the product's sale.

The Warranty does not cover defects resulting from (if applicable):

- User's fault, i.e. product damage caused by unqualified repair work, improper assembly, insufficient insertion of seat post into frame, insufficient tightening of pedals and cranks
- Improper maintenance
- Mechanical damages
- Regular use (e.g. wearing out of rubber and plastic parts, moving mechanisms, joints, wear of brake pads/blocks, chain, tires, cassette/multi wheel etc.)
- Unavoidable event, natural disaster
- Adjustments made by unqualified person
- Improper maintenance, improper placement, damages caused by low or high temperature, water, inappropriate pressure, shocks, intentional changes in design or construction etc.

Warranty Claim Procedure

The Buyer is obliged to check the Goods delivered by the Seller immediately after taking the responsibility for the Goods and its damages, i.e. immediately after its delivery. The Buyer must check the Goods so that he discovers all the defects that can be discovered by such check.

When making a Warranty Claim the Buyer is obliged, on request of the Seller, to prove the purchase and validity of the claim by the Invoice or Bill of Delivery that includes the product's serial number, or eventually by the documents without the serial number. If the Buyer does not prove the validity of the Warranty Claim by these documents, the Seller has the right to reject the Warranty Claim.

If the Buyer gives notice of a defect that is not covered by the Warranty (e.g. in the case that the Warranty Conditions were not fulfilled or in the case of reporting the defect by mistake etc.), the Seller is eligible to require a compensation for all the costs arising from the repair. The cost shall be calculated according to the valid price list of services and transport costs.

If the Seller finds out (by testing) that the product is not damaged, the Warranty Claim is not accepted. The Seller reserves the right to claim a compensation for costs arising from the false Warranty Claim.

In case the Buyer makes a claim about the Goods that is legally covered by the Warranty provided by the Seller, the Seller shall fix the reported defects by means of repair or by the exchange of the damaged part or product for a new one. Based on the agreement of the Buyer, the Seller has the right to exchange the defected Goods for a fully compatible Goods of the same or better technical characteristics. The Seller is entitled to choose the form of the Warranty Claim Procedures described in this paragraph.

The Seller shall settle the Warranty Claim within 30 days after the delivery of the defective Goods, unless a longer period has been agreed upon. The day when the repaired or exchanged Goods is handed over to the Buyer is considered to be the day of the Warranty Claim settlement. When the Seller is not able to settle the Warranty Claim within the agreed period due to the specific nature of the Goods defect, he and the Buyer shall make an agreement about an alternative solution. In case such agreement is not made, the Seller is obliged to provide the Buyer with a financial compensation in the form of a refund.

CZ
SEVEN SPORT s.r.o.

Registered Office: Strakonická 1151/2c, Praha 5, 150 00,
ČR

Headquarters: Dělnická 957, Vítkov, 749 01

Warranty & Service: Čermenská 486, Vítkov 749 01

CRN: 26847264

VAT ID: CZ26847264

Phone: +420 556 300 970

E-mail: eshop@insportline.cz
reklamace@insportline.cz
servis@insportline.cz

Web: www.inSPORTline.cz

SK
inSPORTline s.r.o.

Headquarters, warranty & service center: Električná 6471,
Trenčín 911 01, SK

CRN: 36311723

VAT ID: SK2020177082

Phone: +421(0)326 526 701

E-mail: objednavky@insportline.sk
reklamacie@insportline.sk
servis@insportline.sk

Web: www.inSPORTline.sk

About shipping

